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Operating Instructions

Mode d'emploi

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BOLA Immersion Probes PT 100 LEMO®, Cat.-No.: P 1790-20 / P 1792-20

OPERATING INSTRUCTIONS

1. Supplied with	8
2. Security instructions	8
3. Taking into operation	8
4. Mode of function	9
5. Service / Maintenance	9
6. Disturbances / Repairs	10
7. Waste disposal	10
8. Specifications	10
9. Spare Parts / Accessories	12
10. Customized Temperature Probes	12

BOLA Temperature Probes P 1750-...

BOLA Immersion Probes P 1780-20

BOLA Temperature Probes Lemo® Compact P 1730-...

BOLA Double Temperature Probes Compact P 1740-...

BOLA Temperature Probes Lemo® P 1760-... / P 1762-...

BOLA Immersion Probes Lemo® P 1790-20 / P 1792-20

BOLA Temperature Probes PT 100

Please read the information listed below thoroughly.

We urgently ask you to observe the information about security and use.

BOLA Temperature Probes are PT 100 probes. They are ideal for measuring temperatures in aggressive liquids. They are also available with connectors type LEMO® either plug or socket.

The "LEMO® compact" probe as well as the "Double probe LEMO® compact" have one respective two LEMO® sockets, which are mounted directly on the sensor. The LEMO® connectors are already mounted and can be connected to every suitable measuring instrument.

1. Supplied with

1 Temperature Probe PT 100 with cable length 1.500 mm or Immersion Probe with cable length 4.000mm (both alternatively with or without LEMO® socket / LEMO® plug) respectively Temperature Probe LEMO® Compact (single or double)

1 Operating Instructions


2. Security instructions

The temperature probes are designed to make measurements in aggressive liquids. The operating instructions cannot show all security problems which can arise eventually. It is up to the user to make sure that all safety and health regulations are kept and that the corresponding restrictions are stipulated before use.

- » Please do not take a damaged temperature probe into operation.
- » Please take the response times of the temperature probe into consideration during use. Typical response times T 50 and T 90 (50% respectively 90% of the final measuring value) are to be found under point 8 Specifications on pages 10 and 11.
- » Please do not exceed the operating temperature range of maximum +250 °C and minimum -50 °C.
- » Please verify that the LEMO® connectors are linked correctly to the measuring instrument.

3. Taking into operation

- » The temperature / immersion probes only have to be connected to a suitable measuring instrument. After that, you can start measuring.
- » For the LEMO® compact probes (single and double) please use an adequate cable to connect the probe to the measuring instrument.
- » For temperature and immersion probes without LEMO® connectors, the following connecting plan has to be considered. If you want to mount your own connector, please consider the correct position of the pins.

BOLA Temperature Probes PT 100

4. Mode of function

BOLA temperature probes are sensors based on platinum resistance changes under temperature influence. They are manufactured in conformance to class A. Thus, depending on the real temperature T, the allowed temperature tolerance dT is as follows: the deviation of the measured temperature is equal or less than $\pm 0,15 \text{ }^{\circ}\text{C} + (0,002 \times T)$. The sensor itself is located at the end of a PTFE-encapsulated stainless steel tube (material code 1.4571). This tube provides certain rigidity, but can manually be bent to the requested form.

The seamless PTFE-encapsulation provides an almost universal chemical resistance.

» BOLA Temperature probe PT 100

The probes have tapered tips which reduce the response time considerably. The cable is encapsulated with PFA and connected tightly to the temperature probe. The 4-wire cable has a length of 1.500 mm and is equipped with a LEMO® connector, optionally available as socket size 1 or plug size 1.

» BOLA Immersion probe PT 100


This probe can be immersed completely in aggressive liquids and gases. The 4-wire cable is coated with white PFA. The cable and the probe are resistant to virtually all chemicals. In addition, the total immersion probe can be used for measurement of temperature in autoclaves. They are equipped with a LEMO® connector, optionally available as socket size 1 or plug size 1 and 4-wire cable with a length of 4.000 mm.

» BOLA Probe LEMO® Compact PT 100

The probes have tapered tips which reduce the response time considerably. To connect a cable in any requested length, the LEMO®-socket size 1 has been fixed directly at the end of the probe.

» BOLA Double Probe LEMO® Compact PT 100

This probe is equipped with two independent thermo elements (1) and (2). To connect cables in any requested lengths, the LEMO®-sockets size 1 have been fixed directly at the end of the probe.

5. Service / Maintenance

In case of a strong contamination, the temperature probe should be cleaned after use. Generally, all fluoroplastic materials, PTFE, PFA and FEP, have smooth, non-wetting surfaces and can usually be cleaned without any problems. Please do not use any abrasive scouring agents as they might damage the surfaces. You may use all neutral detergents (pH 7). For a stronger contamination we recommend the use of an alkaline detergent up to pH 12. Apart from that, the probes are maintenance-free.

BOLA Temperature Probes PT 100

6. Disturbances / Repairs

If the probe should not work properly, please check first if the LEMO® connector is linked correctly to the measuring instrument. Also a damage of the connecting cable or of the sensor can lead to disturbances. Should there be a disturbance you cannot handle, please do not hesitate to call us. Of course we will try to help you this way. If we do not find a solution to your problem then, we will ask you to return the cleaned probe in the original packing to us. We will check if a cost-effective repair is possible.

7. Waste disposal

The probes as well as their packing are made of recyclable materials. Please take further information about the materials used from the following specifications.

8. Specifications

» BOLA Temperature Probes PT 100

BOLA Cat. No.	P 1750-...
Materials	Stainless steel tube 1.4571 encapsulated with PTFE, PFA – coated cable
Temperature range	max. +250°C, min. –50°C
Resistance	100 Ohm at 0 °C
Probe-Ø	8 mm
Useful length	100 / 200 / 300 / 500 / 600 mm
Total length, approx	160 / 260 / 360 / 560 / 660 mm
Tip-Ø	5,5 mm
Collar ring	approx. 12 mm
Cable length	1.500 mm
Response time	T 50: 7-12 s; T 90: 14-16 s

» BOLA Temperature Probes PT 100 LEMO® either with socket or plug

BOLA Cat. No.	P 1760-... / P 1762-...
Materials	Stainless steel tube 1.4571 encapsulated with PTFE, PFA – coated cable Chromed brass (LEMO® coupling)
Temperature range	max. +250°C, min. –50°C
Resistance	100 Ohm at 0 °C
Probe-Ø	8 mm
Useful length	100 / 200 / 300 / 500 / 600 mm
Total length, approx.	160 / 260 / 360 / 560 / 660 mm
Tip-Ø	5,5 mm
Collar ring-Ø	approx. 12 mm
Cable length	1.500 mm
Response time	T 50: 7-12 s; T 90: 14-16 s

BOLA Temperature Probes PT 100

» **BOLA Immersion Probes PT 100**

BOLA Cat. No.	P 1780-20
Materials	Stainless steel tube 1.4571 encapsulated with PTFE, PFA – coated cable
Temperature range	max. +250°C, min. -50°C
Resistance	100 Ohm at 0 °C
Probe-Ø	4 mm
Probe length	approx. 50 mm
Cable length	4.000 mm
Response time	T 50: 4-6 s; T 90: 6-8 s

» **BOLA Immersion Probes PT 100 LEMO® either with socket or plug**

BOLA Cat. No.	P 1790-20, P 1792-20
Materials	Stainless steel tube 1.4571 encapsulated with PTFE, PFA – coated cable Chromed brass (LEMO® coupling)
Temperature range	max. +250°C, min. -50°C
Resistance	100 Ohm at 0 °C
Probe-Ø	4 mm
Probe length	approx. 50 mm
Cable length	4.000 mm
Response time	T 50: 4-6 s; T 90: 6-8 s

» **BOLA Temperature Probes PT 100 LEMO® Compact with socket**

BOLA Cat. No.	P 1730-...
Materials	Stainless steel tube 1.4571 encapsulated with PTFE Chromed brass (LEMO® coupling)
Temperature range	max. +250°C, min. -50°C
Resistance	100 Ohm at 0 °C
Probe-Ø	8 mm
Probe length	100/300/500 mm
Total length, approx	170/370/570 mm
Probe tip-Ø	5,5 mm
Collar ring-Ø	ca. 12 mm
Response time	T 50: 7-12 s; T 90: 14-16 s

BOLA Temperature Probes PT 100

» **BOLA Double Temperature Probes PT 100 LEMO® Compact with socket**

BOLA Cat. No.	P 1740-...
Materials	Stainless steel tube 1.4571 encapsulated with PTFE Chromed brass (LEMO® coupling)
Temperature range	max. +250°C, min. -50°C
Resistance	100 Ohm at 0°C
Probe-Ø	8 mm (P 1740-40: 6 mm)
Probe length	300/400/500 mm
Total length, approx	4.000/500/600 mm
Probe tip -Ø	6,0 mm
Collar ring-Ø	27 mm
Response time	T 50: 7-12 s; T 90: 14-16 s

9. Spare Parts / Accessories

» see page 19 and 20

10. Customized Temperature Probes

Should you be in need of a temperature or immersion probe with other dimensions or connectors, simply send us an inquiry with the requested specification. As manufacturers we are the best contact for custom production.

Generally, the following modifications of our standard temperature probes can be offered:

- » Special dimensions of the thermo element (different probe diameter, useful or total length).
- » Different cable lengths..
- » Different kind of LEMO® connectors – a connector of your choice can be mounted.
- » We can coat your temperature probes with a PTFE heat shrinkable tubing so that they have the chemical resistance of PTFE. Because of the thin coating, the probe or thermometer has slower response behaviour. For coating, your probe has to withstand a short-time temperature of minimum +250 °C.

Please do not hesitate to contact us if you have any questions regarding custom manufacture.

» **Registered Trademarks**

LEMO® is a registered trademark by Lemo.

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 **Declaration of Conformity**

EC-Declaration of Conformity in compliance with Machinery Directive 2006/42/EC

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herewith declares that the following products:

Temperature Probes Lemo Compact, product group: P 1730- ...

Double Temperature Probes Lemo Compact, product group: P 1740- ...

Temperature Probes, product group: P 1750- ...

Temperature Probes Lemo, product group: P 1760- ... / P 1762- ...

Total Immersion Probes, product group: P 1780- ...

Total Immersion Probes Lemo, product group: P 1790- ... / P 1792- ...

Temperature Probes PT 1000, product group: P 1950- ...

Temperature Probes PT 1000 Lemo, product group: P 1962- ...

Total Immersion Probes PT 1000, product group: P 1980- ...

Total Immersion Probes PT 1000 Lemo, product group: P 1992- ...

apply to the above mentioned directive including any amendments valid at the time this declaration was signed.

The following directives have been applied:

Directive 2014/30 EU

**Directive relating to electromagnetic compatibility,
General requirements**

The low voltage directive 2006/95/EC has been fulfilled according to 2006/42/EC terms of protection objectives.

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As a manufacturing company we can make products to your specifications. Please send us a drawing or a sketch together with answers to the questions listed below.

Return this form with contact details by fax or e-mail and we will be in touch to discuss details and give you a free, no obligation quote.

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BOHLENDER GmbH